


4. Ion selective electrodes

12 October 2023 09:52


ION SELECTIVE ELECTRODE

The electrode which selectively detects a particular ion (quantitatively measure?)
Eg: glass electrode

Glass Electrode



- pH sensitive electrode
- detects H⁺ ions
- internal reference electrode : Ag | AgCl
- membrane electrode



$$E_b = E_1 - E_2$$

NOTE: Boundary potential
Potential diff. b/w analytical and standard solution

$$E_b = \frac{0.0591}{n} \log \frac{C_1}{C_2}$$

$$E_b = \frac{0.0591}{n} \log C_1 - \frac{0.0591}{n} \log C_2$$

$$E_b = K + \frac{0.0591}{n} \log C_1$$

$$[\because C_1 = [H^+], n = 1]$$

$$E_b = K + 0.0591 \log [H^+]$$

$$[\because pH = -\log[H^+]]$$

$$E_b = K - 0.0591 \text{ pH}$$

CALCULATION OF GLASS ELECTRODE POTENTIAL (E_g)

$$E_g = E_b + E_{\text{reference electrode}} + E_{AP}$$

Asymmetric potential

NOTE: Asymmetric potential
The potential difference b/w

$$E_g = K - 0.0591 \text{ pH} + E_{RE} + E_{AP}$$

$$\therefore E_g^{\circ} = K + E_{RE} + E_{AP},$$

$$E_g = E_g^{\circ} - 0.0591 \text{ pH}$$

Note: Asymmetric potential

The potential difference b/w curvature of glass bulb of analytical soln and standard soln.

APPLICATIONS OF GLASS ELECTRODE

- Industrial
- Agriculture
- Laboratories

ADVANTAGES OF GLASS ELECTRODE

- Can work in both oxidising and reducing environments
- Construction is easy

DISADVANTAGES OF GLASS ELECTRODE

- Sensitive potentiometer since glass bulb is extremely thin
- Alkaline error: glass electrodes using sodium glass bulb (Nall)



Detects only 1-9 pH accurately



10-14 is not detected properly [Alkaline error]

How to reduce alkaline error?


Use lithium glass bulb instead of sodium glass bulb




Accurately detects 1-14 pH but is more expensive

CALCULATION OF pH USING GLASS ELECTRODE

Calomel electrode || Glass electrode





$$E_{\text{cell}} = E_{\text{right}} - E_{\text{left}}$$

$$E_{\text{cell}} = E^\circ_G - 0.0591 \text{ pH} - E_{\text{calomel}}$$

$$E^\circ_G = 0.0591 \text{ pH} + E_{\text{cell}} + E_{\text{calomel}}$$


$$\text{pH} = \frac{E^\circ_G - E_{\text{cell}} - E_{\text{calomel}}}{0.0591}$$

PROBLEMS

Q. A glass electrode is dipped with a solution ($\text{pH} = 4$) offered a potential of 0.284 V at 298 K . With a solution of unknown pH at same temperature, $E = 0.38 \text{ V}$.

$$E_{\text{calomel}} = 0.2422$$

Data



$$\begin{aligned}
 E^\circ_G &= 0.0591(4) + 0.284 + 0.2422 \\
 &= 0.7626
 \end{aligned}$$

$$\begin{aligned}
 \text{pH} &= \frac{E^\circ_G - E_{\text{cell}} - E_{\text{calomel}}}{0.0591} \\
 &= \frac{0.7626 - 0.38 - 0.2422}{0.0591} \\
 &= 2.38
 \end{aligned}$$

Q. Same question: (1) $\text{pH} = 4.066$, $E = 0.209$ at 298 K
 (2) $\text{pH} = ?$, $E = 0.1162$
 $E_{\text{calomel}} = 0.2422$

$$\text{Joh. } E_s^0 = 0.0591 (4.066) + 0.209 + 0.2422 \\ = \underline{\underline{0.6915 \text{ V}}}$$

$$\text{pH} = \frac{E_s^0 - E_{\text{calomel}} - E_{\text{cell}}}{0.0591} \\ = \frac{0.6915 - 0.2422 - 0.1162}{0.0591} \\ = \underline{\underline{5.64}}$$